



The fructophilic yeast for restarting stuck ferments

DESCRIPTION •


UVAFERM 43[™] yeast was selected by the Inter Rhone (France) for its exceptional ability to complete alcoholic fermentations efficiently.

UVAFERM 43[™] was chosen from among 33 different isolates and tested under difficult stuck wine conditions.



BENEFITS & RESULTS

Recent trials that confirmed the UVAFERM 43^{TM} yeast is able to use fructose more readily compared to other wine yeasts. It is the preferred choice for restarting stuck fermentations with high fructose to glucose ratios.




and 65 g/L of glucose)



YSEOTM signifies Yeast Security and Sensory Optimization, a unique Lallemand yeast production process to meet demanding fermentation conditions. While not all yeast benefit from this process, YSEOTM improves the reliability of alcoholic fermentation by improving yeast quality and performance and reduces the risk of organoleptic deviation even under difficult conditions. YSEOTM yeasts are 100% natural and non-GMO.



and 130 g/L of glucose)



- **PROPERTIES*** Saccharomyces cerevisiae Gal. (ex var. bayanus)
 - Alcohol tolerance up to 16% v/v
 - · Good fermentation rate
 - · Competitive factor ("Killer K2") active
 - Low nutritional requirement
 - Low volatile acidity production
 - Low SO₂ and H₂S production
 - · Neutral sensory effect on the finished wine

INSTRUCTIONS FOR OENOLOGICAL USE

A. Rehydration without yeast protector

Dosage rate: 20 to 40 g/hL

- 1. Rehydrate the yeast in 10 times its weight in water (temperature between 35 °C and 40 °C).
- 2. Resuspend the yeast by gently stirring and wait for 20 minutes.
- 3. Mix the rehydrated yeast with a little juice/must, gradually adjusting the yeast suspension temperature to within 5-10 °C of the juice/must temperature.
- 4. Inoculate into the must.

B. Rehydration with a yeast protector

In musts with high alcohol potential (> 13% v/v), with low turbidity (< 80 NTU) or other challenging conditions, the use of one of our GO-FERM™ products (wine yeast protector) during yeast rehydration is recommended. Follow rehydration instructions according to the selected GO-FERM™ product.

Notes:

The total rehydration time should not exceed 45 minutes. It is crucial that a clean container is used to rehydrate the yeast. Rehydration directly in must is generally not advisable. Ensure yeast nutrition is appropriately managed during fermentation.

PACKAGING AND STORAGE

- Available in 500 g and 10 kg
- Store in a cool dry place
- To be used once opened

Distributed by:



1480 Cader Lane, Suite A Petaluma, CA 94954 p. 707-765-6666 | f. 707-765-66 info@scottlab.com | scottlab.com

The information in this document is correct to the best of our knowledge. However, this data sheet should not be considered to be an express guarantee, nor does it have implications as to the sales condition of this product. July 2023.

















^{*}subject to fermentation conditions